

# Sure Vent 2™

## Oxygen Flow Tables

### O<sub>2</sub> Cylinder Autonomy

| Cylinder Size      | D                       | E   | M    | H    |
|--------------------|-------------------------|-----|------|------|
| Capacity in Litres | 387                     | 622 | 3028 | 6905 |
| Flow (LPM)         | Length of Use (minutes) |     |      |      |
| 6                  | 65                      | 100 | 500  | 1150 |
| 8                  | 50                      | 80  | 380  | 860  |
| 10                 | 40                      | 60  | 300  | 690  |
| 12                 | 30                      | 50  | 250  | 570  |
| 15                 | 25                      | 40  | 200  | 460  |
| 20                 | 20                      | 30  | 150  | 340  |
| 25                 | 15                      | 25  | 120  | 270  |
| 30                 | 13                      | 20  | 100  | 230  |
| 35                 | 11                      | 18  | 80   | 190  |
| 40                 | 10                      | 16  | 70   | 170  |

**Entrained Flow table  
for use with SV2™  
optional entainment  
connector**

| 50% Connector     |                      | TOTAL                |
|-------------------|----------------------|----------------------|
| Supply Flow (LPM) | Entrained Flow (LPM) | Delivered Flow (LPM) |
| 6                 | 14                   | 20                   |
| 8                 | 17                   | 25                   |
| 10                | 20                   | 30                   |
| 12                | 23                   | 35                   |
| 15                | 25                   | 40                   |