


ETHICON

ETHILON*

**Monofilament Polyamide
6 or 6/6 Suture**



ETHILON*

Monofilament Polyamide 6 or 6/6 Suture

Monofilament

Extruded to provide a strong uniform strand.

Smooth surface will not support bacterial growth.

Good handling and pliability characteristics permit easy, secure knotting.

Synthetic

For minimal tissue reaction.

Non-Absorbable

Provides prolonged tensile strength retention in tissue.

Colour

Blue or black.

Range










Gauge sizes 11-0 to 2 (U.S.P.). Supplied as needled sutures and ligatures.

Sterilisation

By Irradiation.




These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

ETHILON* Monofilament Polyamide 6 or 6/6 Suture

Needle length and description	Length	Colour	Gauge						
⊙ 1/2 circle ETHIGUARD* Blunt Point			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 <div style="display: flex; justify-content: space-between;"> CTB 40mm </div>	150cm loop	black		W724					
 <div style="display: flex; justify-content: space-between;"> CTXB 48mm </div>	150cm loop	black			W729				
⊙ 1/2 circle Taper Point			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 <div style="display: flex; justify-content: space-between;"> MH-1 31mm </div>	100cm	black	W568		W745				
 <div style="display: flex; justify-content: space-between;"> MH-1 31mm </div>	150cm loop	black			W739				
⊙ 1/2 circle Taper Point (heavier wire diameters)			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 <div style="display: flex; justify-content: space-between;"> CT-1 36mm </div>	150cm loop	black			W743				
 <div style="display: flex; justify-content: space-between;"> CT 40mm </div>	100cm	black		W746	W747				
 <div style="display: flex; justify-content: space-between;"> CT 40mm </div>	150cm loop	black		W740					
 <div style="display: flex; justify-content: space-between;"> CTX 48mm </div>	100cm	black			W749				
 <div style="display: flex; justify-content: space-between;"> CTX heavy 48mm </div>	150cm loop 200cm loop	black black		W741	W748 W768				







ETHILON*



Monofilament Polyamide 6 or 6/6 Suture

Needle length and description	Length	Colour	Gauge				
▼ 1/2 circle TAPERCUT*			2-0 3	0 3.5	1 4	2 5	4 6
 V-39	45mm	200cm loop	black		W757		
▼ 1/2 circle Reverse Cutting			2-0 3	0 3.5	1 4	2 5	4 6
 CP	40mm	100cm	black	W737			
 CPX	48mm	100cm	black		W738		

ETHILON*









Monofilament Polyamide 6 or 6/6 Suture




Needle length and description	Length		Colour	Gauge					
▼ 3/8 circle Reverse Cutting MultiPass* PRIME* Needle				7-0 0.5	6-0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2-0 3.0
 P-1	11mm	45cm	black	MPE1696H	MPE697H				
 P-3	13mm	45cm 45cm	black blue		MPE1698H	MPE698H MPE7800H	MPE699H MPE7801H		
 PS-3	16mm	45cm 45cm 75cm	black blue black		MPE7662H	MPE7930H MPE24012H	MPE7931H		
 PS-2	19mm	45cm 45cm 75cm 75cm	black blue black blue			MPE7804H	MPE1667H MPE1855H MPE24022H MPE7161H	MPE1669H MPE1621H MPE7933H	
 PS-1	24mm	45cm 75cm 75cm	blue black blue				MPE7956H	MPE7805H MPE7563H MPE7954H	MPE7808H
 PS	26mm	45cm 75cm	blue blue					MPE1625H MPE1685H	MPE1626H



▲ 3/8 circle Conventional Cutting MultiPass* PRIME* Needle				7-0 0.5	6-0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2-0 3.0
 PC-3	16mm	45cm 45cm 75cm	black blue blue		MPE1853H	MPE1616H			
 PC-5	19mm	45cm	blue				MPE1856H		

ETHILON*



Monofilament Polyamide 6 or 6/6 Suture


Needle length and description		Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 P-6	8mm	45cm	black	W1600T	W1601T					
 P-1	11mm	45cm 45cm	black blue	W1610T	W1611T	W1612T				
 P-3	13mm	45cm	black	W1851T	W1852T					
 PS-3	16mm	45cm	black	W1614T		W1861T				
 PS-2	19mm	45cm 45cm 75cm	black blue blue		W1618T	W1620T W1619T	W1621T			
 PS-1	24mm	45cm	blue				W1857T			
 PS	26mm	45cm 75cm	blue blue				W1625T W1685T	W1626T		
 PSLX	36mm	75cm	blue				W1627T			





▼ 3/8 circle Reverse Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 FS-2	19mm	45cm	blue			W319				
 FS	26mm	45cm	blue				W320			
 LS-1	45mm	100cm	black					W736		

▲ 3/8 circle Conventional Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 PC-3	16mm	45cm 45cm	black blue	W1615T	W1854T W1616T					
 PC-5	19mm	45cm	blue			W1856T				

ETHILON* Monofilament Polyamide 6 or 6/6 Suture

Needle length and description	Length	Colour	Gauge				
▲ 3/8 circle Conventional Cutting - PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 PC-25 26mm	45cm 45cm	black blue	W1859T	W1858T			
 PCLX 40mm	75cm	blue					W1632T




▲ 3/8 circle Conventional Cutting	Length	Colour	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 CEC135-40 40mm	45cm	blue					W786



▲ 3/8 circle Conventional Cutting (Hand Held)	Length	Colour	2-0 3	0 3.5	1 4	2 5	4 6
 LR-60 CONV 60mm	100cm	blue	W558				
 KP-3 CONV 90mm with 40mm retention tubing	2x50cm	black			W2797		
 KP-3 CONV 90mm	75cm	black			W689		
 KP-3 CONV 90mm	100cm	blue	W795	W796	W797	W798	



ETHILON*


Monofilament Polyamide 6 or 6/6 Suture


Needle length and description				Length	Colour	Gauge				
▼ Straight Cutting						5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5
 KS-55		55mm	100cm	blue					W776	
 KS		60mm	75cm	blue				W788		
 BN		65mm	75cm	blue	W624					


Needle length and description		Length	Colour	Gauge						
 3/8 circle CS ULTIMA* spatula				10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5
 CS140-6	6.5mm	30cm	black	U7061						
 CS140-6	6.5mm	30cm	black	U7003						




 1/2 circle CS ULTIMA* spatula				10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5
 CS160-6	5.5mm	30cm	black	U7000						

▼ 3/8 circle Advanced MICRO-POINT spatula				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG140-6	6.5mm	30cm	black		W1768					
 TG140-6	6.5mm	30cm	black	W1780	W1770	W1769				

▼ 1/2 circle Advanced MICRO-POINT spatula				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG160-6	5.5mm	30cm	black		W1756					















▼ Compound Curve Advanced MICRO-POINT spatula				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG6-C	5mm	30cm	black		W7799					

▼ 3/8 circle MICRO-POINT spatula				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG140-8	6.5mm	30cm	black		W1718	W1717				

▼ 1/2 circle MICRO-POINT spatula				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG160-8	5.5mm	15cm	black		W1721					
 TG175-8	7mm	30cm	black		W1720					
 TG175-8	7mm	30cm	black		W1719					

ETHILON*

Monofilament Polyamide 6 or 6/6 Suture

Needle length and description		Length	Colour	Gauge						
⊙ 3/8 circle Taper Point				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 BV50-3	3.8mm	13cm	black	W2881						
 BV75-3	3.8mm	13cm	black		W2870	W2871				
 BV75-4	5mm	13cm	black		W2850					
 BV100-4	5mm	13cm	black		W2830	W2829				
 BV130-4 VISI-BLACK	5mm	13cm	black			W2913				
 BV130-4	5mm	13cm	black		W2814	W2813	W2812			
 BV130-5 VISI-BLACK	6.5mm	13cm	black				W2908			
 BV130-5	6.5mm	13cm	black		W2810		W2808			
 BV-1	9.3mm	45cm	black				W8170			
⊙ 1/2 circle Taper Point				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 BVH100-3	5mm	13cm	black		W2860					
⊖ 3/8 circle MICRO-POINT				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 V100-3	3.5mm	13cm	black		W2899	W2898				
▾ 3/8 circle MICRO-POINT spatula				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG135-6	4.5mm	13cm	black				W2901			
⊖ 1/4 circle TAPERCUT*				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 V150-4	4.5mm	15cm	black		W1701					
▾ 1/2 circle MICRO-POINT Reverse Cutting				11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 G-7	8mm	45cm	black				W1765			

ETHILON*

Monofilament Polyamide 6 or 6/6 Suture

Ligatures and Non-needed Sutures

Pre-cut Sutures

Gauge		black 2 x 1m lengths 12 sachets per box
1	4	W185
2	5	W186